



1 **Minispir™
Spirometer**
Order Code 911000



2 **Minispir™
Spirometer with Oximeter**
Order Code 911010

Minispir™

Computer based spirometer with oximetry option

✓ **Incomparable Clinical Performance:** Accuracy, Features and EMR Integration/Workflow

Contact your MIR Specialist for Minispir™

FlowMIR™

PATENTED DESIGN

Minispir™ uses the FlowMIR™
cost effective and factory
calibrated disposable
Mouthpiece/Turbine
Calibration Free

FlowMIR™ Turbine prevents
the Minispir™ spirometer
to be exposed to patient
inspired and expired air
Cross-Contamination Free



FlowMIR™ Turbine
Order Code 910004



Contact MIR at 844-4MIRUSA for MIR USA Sales Support & Service, Monday to Friday 8am to 5pm (central time)

Minispir™

- Directly connects to a PC via USB (No **cradle** or **battery** required)
- A great device for Pulmonologists, Respiratory Therapists, Allergists, General Practitioners, Pediatrics, Occupational Medicine, Sports Medicine, and Cardiologists to assess respiratory pathologies

WINSPIROPRO™ SOFTWARE

Each Spirometer includes free WinspiroPRO™ PC Software:

- ✓ Data Management and Interactive Platform that provides incentive routines and enables EMR connectivity
- ✓ includes incentive system to improve patient compliance during Spirometry test

The Minispir™ Spirometer

Includes FVC, VC&IVC, MVV, PRE and POST BD & Bronchial Challenge tests

Wide selection of predicted values

Real time Flow/Volume loop and Volume/Time Curve with PRE/POST comparison

Automatic Spirometry interpretation with alert messages

ATS compliant and supports NHANES III standard

Multi-language interface

Incentive system to improve patient compliance during spirometry test

